

Anti-PEF1 Rabbit antibody



Catalog Number: 168708

Product name

Rabbit IgG

Anti-PEF1 Rabbit antibody

Applications

Specificity

WB, FC, IP

Human, Mouse, Rat

Dilutions

Antibody description

1:1000-1:5000

PEF1 Rabbit monoclonal antibody

1:20

Preparation

1:20-1:50

Antigen: A synthetic peptide of human PEF1

Validations

Formulation

Supplied in 50nM Tris-Glycine(pH 7.4), 0.15M NaCl, 40%Glycerol, 0.01% sodium azide and 0.05% BSA.


Storage

Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze / thaw cycle.

Clonality

Monoclonal

Ig Type



Western blot detection of PEF1 in K562, Rat Brain, C6, 3T3, HeLa cell lysates using PEF1 Rabbit mAb(1:1000 diluted). Predicted band size: 30kDa. Observed band size: 30kDa.