

Anti-Angiotensin II Type 1 Receptor Rabbit antibody



Catalog Number: 168709

Product name

Anti-Angiotensin II Type 1 Receptor Rabbit antibody

Specificity

Human, Mouse, Rat

Antibody description

Angiotensin II Type 1 Receptor Rabbit monoclonal antibody

Preparation

Antigen: A synthetic peptide of human Angiotensin II Type 1 Receptor

Formulation

50nM Tris-Glycine(pH 7.4), 0.15M NaCl, 40%Glycerol, 0.01% sodium azide and 0.05% BSA.

Storage

Store at 4°C short term. Store at -20°C long term. Avoid freeze / thaw cycle.

Clonality

Monoclonal

Ig Type

Rabbit IgG


Applications

WB

Dilutions

WB: 1/1000-1/5000

Validations



Western blot detection of Angiotensin II Type 1 Receptor in HepG2 cell lysates using Angiotensin II Type 1 Receptor Rabbit mAb(1:1000 diluted).Predicted band size:41kDa.Observed band size:41kDa.