

Anti-Legumain Rabbit antibody



Catalog Number: 168712

Product name

Anti-Legumain Rabbit antibody

Specificity

Human, Mouse, Rat

Antibody description

Legumain Rabbit monoclonal antibody

Preparation

Antigen: Recombinant protein of human Legumain

Formulation

50nM Tris-Glycine(pH 7.4), 0.15M NaCl, 40%Glycerol, 0.01% sodium azide and 0.05% BSA.

Storage

Store at 4°C short term. Store at -20°C long term. Avoid freeze / thaw cycle.

Clonality

Monoclonal

Ig Type

Rabbit IgG


Applications

WB

Dilutions

WB: 1/2000-1/10000

Validations



Western blot detection of Legumain in Jurkat, HeLa cell lysates using Legumain Rabbit mAb(1:1000 diluted). Predicted band size: 49kDa. Observed band size: 56, 46, 37kDa.