

Anti-DBF4 Rabbit antibody



Catalog Number: 168713

Product name

Anti-DBF4 Rabbit antibody

Specificity

Human

Antibody description

DBF4 Rabbit monoclonal antibody

Preparation

Antigen: A synthetic peptide of human DBF4

Formulation

50nM Tris-Glycine(pH 7.4), 0.15M NaCl, 40%Glycerol, 0.01% sodium azide and 0.05% BSA.

Storage

Store at 4°C short term. Store at -20°C long term. Avoid freeze / thaw cycle.

Clonality

Monoclonal

Ig Type

Rabbit IgG


Applications

WB

Dilutions

WB: 1/2000-1/10000

Validations



Western blot detection of DBF4 in HeLa, Jurkat cell lysates using DBF4 Rabbit mAb(1:500 diluted). Predicted band size: 77kDa. Observed band size: 77kDa.