

# Anti-TRIM21 Rabbit antibody



Catalog Number: 168716

## Product name

Anti-TRIM21 Rabbit antibody

## Specificity

Human, Mouse, Rat

## Antibody description

TRIM21 Rabbit monoclonal antibody

## Preparation

Antigen: Recombinant protein of human TRIM21

## Formulation

Supplied in 50nM Tris-Glycine(pH 7.4), 0.15M NaCl, 40%Glycerol, 0.01% sodium azide and 0.05% BSA.

## Storage

Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze / thaw cycle.

## Clonality

Monoclonal

## Ig Type

Rabbit IgG


## Applications

WB

## Dilutions

1:1000

## Validations



Western blot detection of TRIM21 in K562, HeLa cell lysates using TRIM21 Rabbit mAb(1:1000 diluted). Predicted band size: 54kDa. Observed band size: 50kDa.