

# Anti-14-3-3 Rabbit antibody



Catalog Number: 168725

## Product name

Anti-14-3-3 Rabbit antibody

## Specificity

Human, Mouse, Rat

## Antibody description

14-3-3 Rabbit monoclonal antibody

## Preparation

Antigen: A synthetic peptide of human 14-3-3

## Formulation

50nM Tris-Glycine(pH 7.4), 0.15M NaCl,  
40%Glycerol, 0.01% sodium azide and 0.05% BSA.

## Storage

Store at 4°C short term. Store at -20°C long term.  
Avoid freeze / thaw cycle.

## Clonality

Monoclonal

## Ig Type

Rabbit IgG


## Applications

WB

## Dilutions

WB: 1/2000

## Validations



Western blot detection of 41701 in Jurkat,Rat Brain,C6,CHO-K1 cell lysates using 14-3-3 Rabbit mAb(1:500 diluted).Predicted band size:28kDa.Observed band size:27-29kDa.