

Anti-Caveolin-1 Rabbit antibody



Catalog Number: 168735

Product name

Anti-Caveolin-1 Rabbit antibody

Specificity

Human, Mouse

Antibody description

Caveolin-1 Rabbit monoclonal antibody

Preparation

Antigen: Recombinant protein of human Caveolin-1

Formulation

Supplied in 50nM Tris-Glycine(pH 7.4), 0.15M NaCl, 40%Glycerol, 0.01% sodium azide and 0.05% BSA.

Storage

Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze / thaw cycle.

Clonality

Monoclonal

Ig Type

Rabbit IgG

Applications

WB, IHC-P, ICC/IF, FC, IP

Dilutions

1:2000


1:20-1:500

1:50

1:20

1:20

Validations



Western blot detection of Caveolin-1 in HeLa cell lysates using Caveolin-1 Rabbit mAb(1:1000 diluted).Predicted band size:20kDa.Observed band size:20kDa.