

# Anti-StAR Rabbit antibody



Catalog Number: 168738

## Product name

Anti-StAR Rabbit antibody

## Specificity

Human

## Antibody description

StAR Rabbit monoclonal antibody

## Preparation

Antigen: A synthetic peptide of human StAR

## Formulation

Supplied in 50nM Tris-Glycine(pH 7.4), 0.15M NaCl, 40%Glycerol, 0.01% sodium azide and 0.05% BSA.

## Storage

Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze / thaw cycle.

## Clonality

Monoclonal

## Ig Type

Rabbit IgG

## Applications


WB, IHC-P

## Dilutions

1:2000

1:50

## Validations



Western blot detection of StAR in K562,3T3,HeLa cell lysates using StAR Rabbit mAb(1:1000 diluted).Predicted band size:32kDa.Observed band size:32kDa.