

Anti-Xanthine Oxidase Rabbit antibody



Catalog Number: 168761

Product name

Anti-Xanthine Oxidase Rabbit antibody

Specificity

Human, Mouse, Rat

Antibody description

Xanthine Oxidase Rabbit monoclonal antibody

Preparation

Antigen: A synthetic peptide of human Xanthine Oxidase

Formulation

Supplied in 50nM Tris-Glycine(pH 7.4), 0.15M NaCl, 40%Glycerol, 0.01% sodium azide and 0.05% BSA.

Storage

Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze / thaw cycle.

Clonality

Monoclonal

Ig Type

Rabbit IgG


Applications

WB

Dilutions

WB: 1/1000

Validations



Western blot detection of Xanthine Oxidase in K562 cell lysates using Xanthine Oxidase Rabbit mAb(1:1000 diluted).Predicted band size:146kDa.Observed band size:146kDa.