

# Anti-NOLC1 Rabbit antibody



Catalog Number: 168767

## Product name

Anti-NOLC1 Rabbit antibody

## Specificity

Human, Mouse

## Antibody description

NOLC1 Rabbit monoclonal antibody

## Preparation

Antigen: A synthetic peptide of human NOLC1

## Formulation

50nM Tris-Glycine(pH 7.4), 0.15M NaCl, 40%Glycerol, 0.01% sodium azide and 0.05% BSA.

## Storage

Store at 4°C short term. Store at -20°C long term. Avoid freeze / thaw cycle.

## Clonality

Monoclonal

## Ig Type

Rabbit IgG

## Applications

WB, IHC-P, ICC/IF, FC

## Dilutions


WB: 1/2000-1/10000

IHC: 1/20

ICC/IF: 1/100

FC: 1/50

## Validations



Western blot detection of NOLC1 in Jurkat cell lysates using NOLC1 Rabbit mAb(1:500 diluted).Predicted band size:74kDa.Observed band size:110kDa.