

# Anti-Phospho-eIF4EBP1 (Thr46) Rabbit antibody



Catalog Number: 168773

## Product name

Rabbit IgG

Anti-Phospho-eIF4EBP1 (Thr46) Rabbit antibody

## Applications

## Specificity

WB, IHC-P, ICC/IF, IP

Human, Mouse

## Dilutions

## Antibody description

1:2000

Phospho-eIF4EBP1 (Thr46) Rabbit monoclonal antibody

1:50

## Preparation

1:50

1:20

Antigen: A synthetic phosphopeptide corresponding to residues surrounding Thr45 of human eIF4EBP1

## Validations

## Formulation

Supplied in 50nM Tris-Glycine(pH 7.4), 0.15M NaCl, 40%Glycerol, 0.01% sodium azide and 0.05% BSA.


## Storage

Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze / thaw cycle.

## Clonality

Monoclonal

## Ig Type



Western blot detection of Phospho-eIF4EBP1 (Thr46) in K562, C6, 3T3, HeLa cell lysates using Phospho-eIF4EBP1 (Thr46) Rabbit mAb(1:1000 diluted). Predicted band size: 13kDa. Observed band size: 15-20kDa.