

Anti-Cytokeratin 18 Rabbit antibody



Catalog Number: 168790

Product name

Rabbit IgG

Anti-Cytokeratin 18 Rabbit antibody

Applications

Specificity

WB, IHC-P, ICC/IF, FC

Human

Dilutions

Antibody description

WB: 1/2000

Cytokeratin 18 Rabbit monoclonal antibody

IHC: 1/20

Preparation

ICC/IF: 1/20-1/100


Antigen: Recombinant protein of human Cytokeratin 18

FC: 1/20

Formulation

50nM Tris-Glycine(pH 7.4), 0.15M NaCl, 40%Glycerol, 0.01% sodium azide and 0.05% BSA.

Validations



Storage

Store at 4°C short term. Store at -20°C long term. Avoid freeze / thaw cycle.

Western blot detection of Cytokeratin 18 in HeLa cell lysates using Cytokeratin 18 Rabbit mAb(1:500 diluted).Predicted band size:48kDa.Observed band size:48kDa.

Clonality

Monoclonal

Ig Type