

Anti-HDAC2 Rabbit antibody



Catalog Number: 168794

Product name 1:1000

Anti-HDAC2 Rabbit antibody 1:500

Specificity 1:200

Human, Mouse, Rat 1:50

Antibody description 1:20

HDAC2 Rabbit monoclonal antibody

Preparation

Antigen: Recombinant protein of human HDAC2

Formulation

Supplied in 50nM Tris-Glycine(pH 7.4), 0.15M NaCl, 40%Glycerol, 0.01% sodium azide and 0.05% BSA.

Storage

Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze / thaw cycle.

Clonality

Monoclonal

Ig Type


Rabbit IgG

Applications


WB, IHC-P, ICC/IF, FC, IP

Dilutions

Validations



Western blot detection of HDAC2 in K562, C6, 3T3 cell lysates using HDAC2 Rabbit mAb(1:1000 diluted). Predicted band size: 55kDa. Observed band size: 60kDa.



Immunocytochemistry of HDAC2(green) in HeLa cells using HDAC2 Rabbit mAb at dilution 1/200, and DAPI(blue)