

# Anti-HDAC1 Rabbit antibody



Catalog Number: 168801

## Product name

Anti-HDAC1 Rabbit antibody

## Specificity

Human, Mouse, Rat

## Antibody description

HDAC1 Rabbit monoclonal antibody

## Preparation

Antigen: A synthetic peptide of human HDAC1

## Formulation

Supplied in 50nM Tris-Glycine(pH 7.4), 0.15M NaCl, 40%Glycerol, 0.01% sodium azide and 0.05% BSA.

## Storage

Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze / thaw cycle.

## Clonality

Monoclonal

## Ig Type

Rabbit IgG

## Applications

WB, ICC/IF, FC


## Dilutions

WB: 1/2000

ICC/IF: 1/50

FC: 1/20

## Validations



Western blot detection of HDAC1 in K562 cell lysates using HDAC1 Rabbit mAb(1:1000 diluted). Predicted band size:55kDa.Observed band size:62kDa.