

# Anti-Phospho-ERK1 (Thr202)/ERK2 (Thr185) Rabbit antibody



Catalog Number: 168811

## Product name

Anti-Phospho-ERK1 (Thr202)/ERK2 (Thr185)  
Rabbit antibody

## Specificity

Human, Mouse, Rat

## Antibody description

Phospho-ERK1 (Thr202)/ERK2 (Thr185) Rabbit  
monoclonal antibody

## Preparation

Antigen: A synthetic phosphopeptide  
corresponding to residues surrounding Thr185 of  
human ERK2

## Formulation

Supplied in 50nM Tris-Glycine(pH 7.4), 0.15M  
NaCl, 40%Glycerol, 0.01% sodium azide and  
0.05% BSA.

## Storage

Store at 4°C short term. Aliquot and store at -20°C  
long term. Avoid freeze / thaw cycle.

## Clonality

Monoclonal

## Ig Type

Rabbit IgG

## Applications

WB, IHC-P, ICC/IF, IP

## Dilutions


1:1000

1:20

1:20

1:20

## Validations



Western blot detection of Phospho-p44/42 MAPK  
(ERK1/2) (Thr202) in C6,CHO-K1 cell lysates using  
Phospho-p44/42 MAPK (ERK1/2) (Thr202) Rabbit  
mAb(1:1000 diluted).Predicted band  
size:41kDa.Observed band size:44,42kDa.