

Anti-Phospho-ERK1 (Thr202)/ERK2 (Thr185) Rabbit antibody



Catalog Number: 168811

Product name

Anti-Phospho-ERK1 (Thr202)/ERK2 (Thr185)
Rabbit antibody

Ig Type

Rabbit IgG

Specificity

Human, Mouse, Rat

Applications

WB, IHC-P, ICC/IF, IP

Antibody description

Phospho-ERK1 (Thr202)/ERK2 (Thr185) Rabbit
monoclonal antibody

Dilutions

1:1000

1:20

1:20

1:20

Preparation

Antigen: A synthetic phosphopeptide
corresponding to residues surrounding Thr185 of
human ERK2

Formulation

Supplied in 50nM Tris-Glycine(pH 7.4), 0.15M
NaCl, 40%Glycerol, 0.01% sodium azide and
0.05% BSA.


Storage

Store at 4°C short term. Aliquot and store at -20°C
long term. Avoid freeze / thaw cycle.

Clonality

Monoclonal

Validations



Western blot detection of Phospho-p44/42 MAPK
(ERK1/2) (Thr202) in C6,CHO-K1 cell lysates using
Phospho-p44/42 MAPK (ERK1/2) (Thr202) Rabbit
mAb(1:1000 diluted).Predicted band
size:41kDa.Observed band size:44,42kDa.