

# Anti-JNK1 Rabbit antibody

Catalog Number: 168812

## Product name

Anti-JNK1 Rabbit antibody

## Specificity

Human, Mouse, Rat

## Antibody description

JNK1 Rabbit monoclonal antibody

## Preparation

Antigen: Recombinant protein of human JNK1

## Formulation

Supplied in 50nM Tris-Glycine(pH 7.4), 0.15M NaCl, 40%Glycerol, 0.01% sodium azide and 0.05% BSA.

## Storage

Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze / thaw cycle.

## Clonality

Monoclonal

## Ig Type

Rabbit IgG


## Applications

WB

## Dilutions

1:1000

## Validations



Western blot detection of JNK1 in Jurkat,Rat Brain,C6 cell lysates using JNK1 Rabbit mAb(1:1000 diluted).Predicted band size:48kDa.Observed band size:46kDa.