

Anti-MYO19 Rabbit antibody



Catalog Number: 168814

Product name

Anti-MYO19 Rabbit antibody

Specificity

Human

Antibody description

MYO19 Rabbit monoclonal antibody

Preparation

Antigen: Recombinant protein of human MYO19

Formulation

Supplied in 50nM Tris-Glycine(pH 7.4), 0.15M NaCl, 40%Glycerol, 0.01% sodium azide and 0.05% BSA.

Storage

Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze / thaw cycle.

Clonality

Monoclonal

Ig Type

Rabbit IgG

Applications


WB, ICC/IF

Dilutions

1:1000

1:20

Validations



Western blot detection of MYO19 in K562,3T3,HeLa cell lysates using MYO19 Rabbit mAb(1:1000 diluted).Predicted band size:109kDa.Observed band size:109kDa.