

# Anti-Ikaros Rabbit antibody



Catalog Number: 168839

## Product name

Anti-Ikaros Rabbit antibody

## Specificity

Human, Mouse, Rat

## Antibody description

Ikaros Rabbit monoclonal antibody

## Preparation

Antigen: A synthetic peptide of human Ikaros

## Formulation

50nM Tris-Glycine(pH 7.4), 0.15M NaCl, 40%Glycerol, 0.01% sodium azide and 0.05% BSA.

## Storage

Store at 4°C short term. Store at -20°C long term. Avoid freeze / thaw cycle.

## Clonality

Monoclonal

## Ig Type

Rabbit IgG

## Applications

WB, IHC-P, FC


## Dilutions

WB: 1/2000

IHC: 1/200

FC: 1/20

## Validations



Western blot detection of Ikaros in CHO-K1, Jurkat cell lysates using Ikaros Rabbit mAb(1:500 diluted). Predicted band size: 58kDa. Observed band size: 50-70kDa.