

# Anti-Protein Kinase A regulatory subunit I alpha/PRKAR1A Rabbit antibody



Catalog Number: 168842

## Product name

Rabbit IgG

Anti-Protein Kinase A regulatory subunit I alpha/PRKAR1A Rabbit antibody

## Applications

WB, IHC-P, ICC/IF

## Specificity

Human, Mouse, Rat

## Dilutions

WB: 1/1000-1/5000

## Antibody description

Protein Kinase A regulatory subunit I alpha/PRKAR1A Rabbit monoclonal antibody

IHC: 1/20

ICC/IF: 1/20


## Preparation

Antigen: A synthetic peptide of human Protein Kinase A regulatory subunit I alpha

## Validations

## Formulation

50nM Tris-Glycine(pH 7.4), 0.15M NaCl, 40%Glycerol, 0.01% sodium azide and 0.05% BSA.



## Storage

Store at 4°C short term. Store at -20°C long term. Avoid freeze / thaw cycle.

Western blot detection of Protein Kinase A regulatory subunit I alpha/PRKAR1A in C6,Rat Brain,Jurkat cell lysates using Protein Kinase A regulatory subunit I alpha/PRKAR1A Rabbit mAb(1:500 diluted).Predicted band size:43kDa.Observed band size:43kDa.

## Clonality

Monoclonal

## Ig Type