

# Anti-R-Cadherin/Cadherin 4 Rabbit antibody



Catalog Number: 168852

## Product name

Anti-R-Cadherin/Cadherin 4 Rabbit antibody

## Ig Type

Rabbit IgG

## Specificity

Human, Mouse, Rat

## Applications

WB, IP

## Antibody description

R-Cadherin/Cadherin 4 Rabbit monoclonal antibody

## Dilutions

1:2000

1:20

## Preparation

Antigen: A synthetic peptide of human R Cadherin

## Validations

## Formulation


Supplied in 50nM Tris-Glycine(pH 7.4), 0.15M NaCl, 40%Glycerol, 0.01% sodium azide and 0.05% BSA.

## Storage

Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze / thaw cycle.

## Clonality

Monoclonal



Western blot detection of R-Cadherin/Cadherin 4 in Rat Brain, C6 cell lysates using R-Cadherin/Cadherin 4 Rabbit mAb(1:1000 diluted). Predicted band size: 100kDa. Observed band size: 130kDa.