

# Anti-MLH1 Rabbit antibody



Catalog Number: 168860

## General Information

### Product name:

Anti-MLH1 Rabbit antibody

### Immunogen:

### Applications:

WB

### Specificity:

Human, Mouse, Rat

### Antibody description:

MLH1 Rabbit monoclonal antibody

## Properties

### Preparation:

Antigen: A synthetic peptide of human MLH1

### Formulation:

Supplied in 50nM Tris-Glycine(pH 7.4), 0.15M NaCl, 40%Glycerol, 0.01% sodium azide and 0.05% BSA.

### Storage:

Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze / thaw cycle.

### Clonality:


Monoclonal antibody

Ig Type:Rabbit IgG

## Applications

1:1000-1:5000

## Validations



Western blot detection of MLH1 in Rat Brain lysates using MLH1 Rabbit mAb(1:1000 diluted).Predicted band size:85kDa.Observed band size:85kDa.