

Anti-hnRNP G-T Rabbit antibody



Catalog Number: 168875

General Information

Product name:

Anti-hnRNP G-T Rabbit antibody

Immunogen:

Applications:

WB

Specificity:

Human, Mouse, Rat

Antibody description:

hnRNP G-T Rabbit monoclonal antibody

Properties

Preparation:

Antigen: Recombinant protein of human hnRNP G-T

Formulation:

Supplied in 50nM Tris-Glycine(pH 7.4), 0.15M NaCl, 40%Glycerol, 0.01% sodium azide and 0.05% BSA.

Storage:

Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze / thaw cycle.

Clonality:


Monoclonal antibody

Ig Type:Rabbit IgG

Applications

WB: 1/1000

Validations



Western blot detection of hnRNP G-T in K562,Rat Brain,C6,3T3,HeLa cell lysates using hnRNP G-T Rabbit mAb(1:1000 diluted).Predicted band size:43kDa.Observed band size:43kDa.