

Anti-QK1 Rabbit antibody



Catalog Number: 168884

Product name

Anti-QK1 Rabbit antibody

Specificity

Human, Mouse, Rat

Antibody description

QK1 Rabbit monoclonal antibody

Preparation

Antigen: A synthetic peptide of human QK1

Formulation

50nM Tris-Glycine(pH 7.4), 0.15M NaCl, 40%Glycerol, 0.01% sodium azide and 0.05% BSA.

Storage

Store at 4°C short term. Store at -20°C long term. Avoid freeze / thaw cycle.

Clonality

Monoclonal

Ig Type

Rabbit IgG

Applications

WB, ICC/IF, IP


Dilutions

WB: 1/1000-1/5000

ICC/IF: 1/20-1/50

IP: 1/20-1/50

Validations



Western blot detection of QK1 in K562 cell lysates using QK1 Rabbit mAb(1:1000 diluted). Predicted band size:38kDa.Observed band size:38kDa.