

Anti-Glutaminase C Rabbit antibody



Catalog Number: 168897

General Information

Product name:

Anti-Glutaminase C Rabbit antibody

Immunogen:

Applications:

WB, ICC/IF, FC, IP

Specificity:

Human, Mouse, Rat

Antibody description:

Glutaminase C Rabbit monoclonal antibody

Properties

Preparation:

Antigen: Recombinant protein of human
Glutaminase C

Formulation:

Supplied in 50nM Tris-Glycine(pH 7.4), 0.15M
NaCl, 40%Glycerol, 0.01% sodium azide and
0.05% BSA.

Storage:

Store at 4°C short term. Aliquot and store at -20°C
long term. Avoid freeze / thaw cycle.

Clonality:

Monoclonal antibody

Ig Type:Rabbit IgG

Applications


1:1000

1:50

1:50

1:20

Validations



Western blot detection of Glutaminase in
K562,3T3,HeLa cell lysates using Glutaminase
Rabbit mAb(1:1000 diluted).Predicted band
size:74kDa.Observed band size:68kDa.