

Anti-CD99 Rabbit antibody



Catalog Number: 168931

General Information

Product name:

Anti-CD99 Rabbit antibody

Immunogen:

Applications:

WB, IHC-P, ICC/IF, FC, IP

Specificity:

Human

Antibody description:

CD99 Rabbit monoclonal antibody

Properties

Preparation:

Antigen: A synthetic peptide of human CD99

Formulation:

Supplied in 50nM Tris-Glycine(pH 7.4), 0.15M NaCl, 40%Glycerol, 0.01% sodium azide and 0.05% BSA.

Storage:

Store at 4°C short term. Aliquot and store at -20°C

long term. Avoid freeze / thaw cycle.

Clonality:

Monoclonal antibody

Ig Type:

Applications

1:1000-1:5000


1:200-1:500

1:50-1:200

1:20

1:20

Validations



Western blot detection of CD99 in 3T3 cell lysates using CD99 Rabbit mAb(1:1000 diluted). Predicted band size:19kDa.Observed band size:30kDa.