

Anti-CD147 Rabbit antibody



Catalog Number: 168976

Product name

Anti-CD147 Rabbit antibody

Specificity

Mouse

Antibody description

CD147 Rabbit monoclonal antibody

Preparation

Antigen: Recombinant protein of mouse CD147

Formulation

Supplied in 50nM Tris-Glycine(pH 7.4), 0.15M NaCl, 40%Glycerol, 0.01% sodium azide and 0.05% BSA.

Storage

Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze / thaw cycle.

Clonality

Monoclonal

Ig Type

Rabbit IgG

Applications

WB, ICC/IF, FC, IP

Dilutions


WB: 1/2000

ICC/IF: 1/50

FC: 1/50

IP: 1/20

Validations



Western blot detection of CD147 in Rat Brain, C6, 3T3 cell lysates using CD147 Rabbit mAb(1:1000 diluted). Predicted band size: 43kDa. Observed band size: 38-58kDa.