

Anti-ENO1 Rabbit antibody

Catalog Number: 168978

Product name

Anti-ENO1 Rabbit antibody

Dilutions

1:1000

Specificity

1:500

Human, Mouse, Rat

1:100

Antibody description

1:100

ENO1 Rabbit monoclonal antibody

1:20

Validations

Preparation

Antigen: A synthetic peptide of human ENO1

Formulation

Supplied in 50nM Tris-Glycine(pH 7.4), 0.15M NaCl, 40%Glycerol, 0.01% sodium azide and 0.05% BSA.

Storage

Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze / thaw cycle.

Clonality


Monoclonal

Ig Type


Rabbit IgG

Applications

WB, IHC-P, ICC/IF, FC, IP



Immunocytochemistry of ENO1(green) in HeLa cells using ENO1 Rabbit mAb at dilution 1/50, and



Western blot detection of ENO1 in Jurkat,Rat Brain,C6,HeLa cell lysates using ENO1 Rabbit mAb(1:1000 diluted).Predicted band size:47kDa.Observed band size:47kDa.