

Anti-Tau Rabbit antibody



Catalog Number: 168980

General Information

Product name:

Anti-Tau Rabbit antibody

Immunogen:

Applications:

WB, IHC-P, IP

Specificity:

Human, Mouse, Rat

Antibody description:

Tau Rabbit monoclonal antibody

Properties

Preparation:

Antigen: A synthetic peptide of human Tau

Formulation:

Supplied in 50nM Tris-Glycine(pH 7.4), 0.15M NaCl, 40%Glycerol, 0.01% sodium azide and 0.05% BSA.

Storage:

Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze / thaw cycle.

Clonality:

Monoclonal antibody

Ig Type:Rabbit IgG


Applications

1:1000

1:20-1:500

1:20

Validations



Western blot detection of Tau in Rat Brain cell using Tau Rabbit mAb(1:1000 diluted).Predicted band size:79kDa.Observed band size:50-80kDa.