

# Anti-eIF1A Rabbit antibody



Catalog Number: 168992

## General Information

### Product name:

Anti-eIF1A Rabbit antibody

### Immunogen:

### Applications:

WB, IHC-P, ICC/IF, IP

### Specificity:

Human, Mouse, Rat

### Antibody description:

eIF1A Rabbit monoclonal antibody

## Properties

### Preparation:

Antigen: A synthetic peptide of human eIF1A

### Formulation:

50nM Tris-Glycine(pH 7.4), 0.15M NaCl, 40%Glycerol, 0.01% sodium azide and 0.05% BSA.

### Storage:

Store at 4°C short term. Store at -20°C long term. Avoid freeze / thaw cycle.

### Clonality:

Monoclonal antibody

**Ig Type:**Rabbit IgG

## Applications


WB: 1/1000

IHC: 1/20

ICC/IF: 1/20

IP: 1/20

## Validations



Western blot detection of eIF1A in LnCap, K562, U251MG, HepG2, MCF7, HeLa cell lysates using eIF1A Rabbit mAb(1:500 diluted). Predicted band size: 17kDa. Observed band size: 17kDa.