

# Anti-UFM1 Rabbit antibody

Catalog Number: 168999

## Product name

Anti-UFM1 Rabbit antibody

## Specificity

Human, Mouse, Rat

## Antibody description

UFM1 Rabbit monoclonal antibody

## Preparation

Antigen: A synthetic peptide of human UFM1

## Formulation

Supplied in 50nM Tris-Glycine(pH 7.4), 0.15M NaCl, 40%Glycerol, 0.01% sodium azide and 0.05% BSA.

## Storage

Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze / thaw cycle.

## Clonality

Monoclonal

## Ig Type

Rabbit IgG

## Applications

WB, IHC-P, IP


## Dilutions

1:1000


1:50-1:200

1:20

## Validations



Western blot detection of UFM1 in Jurkat cell lysates using UFM1 Rabbit mAb(1:1000 diluted). Predicted band size:9kDa. Observed band



size:9kDa.

Immunohistochemistry of UFM1 in paraffin-embedded human lung cancer tissue using UFM1 Rabbit mAb at dilution 1/100