

Anti-Cathepsin B Rabbit antibody



Catalog Number: 169010

Product name

Anti-Cathepsin B Rabbit antibody

Specificity

Mouse, Rat

Antibody description

Cathepsin B Rabbit monoclonal antibody

Preparation

Antigen: Recombinant protein of mouse Cathepsin B

Formulation

Supplied in 50nM Tris-Glycine(pH 7.4), 0.15M NaCl, 40%Glycerol, 0.01% sodium azide and 0.05% BSA.

Storage

Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze / thaw cycle.

Clonality

Monoclonal

Ig Type

Rabbit IgG

Applications

WB, ICC/IF, FC


Dilutions

WB: 1/1000

ICC/IF: 1/50

FC: 1/100

Validations



Western blot detection of Cathepsin B in Rat Brain, C6, 3T3 cell lysates using Cathepsin B Rabbit mAb(1:1000 diluted). Predicted band size: 37kDa. Observed band size: 44, 27, 24kDa.