

# Anti-TAGAP Rabbit antibody



Catalog Number: 169023

## Product name

Anti-TAGAP Rabbit antibody

## Specificity

Human

## Antibody description

TAGAP Rabbit monoclonal antibody

## Preparation

Antigen: Recombinant protein of human TAGAP

## Formulation

50nM Tris-Glycine(pH 7.4), 0.15M NaCl,  
40%Glycerol, 0.01% sodium azide and 0.05% BSA.

## Storage

Store at 4°C short term. Store at -20°C long term.  
Avoid freeze / thaw cycle.

## Clonality

Monoclonal

## Ig Type

Rabbit IgG

## Applications

WB, FC, IP


## Dilutions

WB: 1/1000

FC: 1/20

IP: 1/20

## Validations



Western blot detection of TAGAP in Rat Brain,C6,CHO-K1,HeLa cell lysates using TAGAP Rabbit mAb(1:1000 diluted).Predicted band size:81kDa.Observed band size:81kDa.