

Anti-Phospho-p44/42 MAPK (ERK1/2) (Thr202/Tyr204) Rabbit antibody



Catalog Number: 169073

General Information

Product name:

Anti-Phospho-p44/42 MAPK (ERK1/2)
(Thr202/Tyr204) Rabbit antibody

Immunogen:

Applications:

WB, IP

Specificity:

Human, Mouse, Rat

Antibody description:

Phospho-p44/42 MAPK (ERK1/2) (Thr202/Tyr204)
Rabbit polyclonal antibody

Properties

Preparation:

Antigen: A synthesized peptide derived from
human Phospho-ERK1 (Y204) + ERK2 (Y187)

Formulation:

Rabbit IgG in phosphate buffered saline , pH 7.4,
150mM NaCl, 0.02% sodium azide and 50%
glycerol.

Storage:

Store at 4°C short term. Store at -20°C long term.
Avoid freeze / thaw cycle.

Clonality:

Polyclonal antibody


Ig Type:Rabbit IgG

Applications

WB: 1/500-1/2000

IP: 1/50

Validations



Western blot analysis of Phospho-Erk1
(T202/Y204) + Erk2 (T185/Y187) expression in
A431 cell lysate treated with EGF.