

Anti-GAPDH Mouse mAb (HRP Conjugate)



Catalog Number: 169100

General Information

Product name:

Anti-GAPDH Mouse mAb (HRP Conjugate)

Immunogen:

Applications:

WB, ELISA

Specificity:

Human, Mouse, Rat

Antibody description:

GAPDH Mouse mAb (HRP Conjugate)

Properties

Preparation:

Antigen: A synthetic peptide corresponding to residues near the amino terminus of human GAPDH protein.

Formulation:

In PBS(pH 7.4) containing with 50% glycerol.

Storage:

Store at 4°C short term. Store at -20°C long term. Avoid freeze / thaw cycle.

Clonality:


Monoclonal

Ig Type: IgG1

Applications

WB: 1/500—1/5000

Validations



Western blot detection of GAPDH in Rat brain, CHO-K1 and HUVEC cell lysates using GAPDH Mouse mAb (HRP Conjugate) (1:5000 diluted). Predicted band size: 36kDa. Observed band size: 36kDa.