

# Anti-IL-1 beta Rabbit antibody



Catalog Number: 169102

## Product name

Anti-IL-1 beta Rabbit antibody

## Specificity

Human, Mouse, Rat

## Antibody description

IL-1 beta Rabbit polyclonal antibody

## Preparation

Antigen: Recombinant fusion protein containing a sequence corresponding to amino acids 124-223 of human IL1B (NP\_000567.1).CTLRDSQQKSLVMSGPYELKALHLQGQDMEQQVVFMSFVQGEESNDKIPVALGLKEKNLYLSCVLKDDKPTLQLESVDPKNYPKKKMEKRFVFNKIEIN

## Formulation

PBS with 0.02% sodium azide,50% glycerol,pH7.3.

## Storage

Store at -20°C . Avoid freeze / thaw cycles.

## Clonality

Polyclonal

## Ig Type

Rabbit IgG


## Applications

WB


## Dilutions

WB: 1/500 - 1/2000

## Validations



Western blot - IL1B Polyclonal Antibody. Western blot analysis of extracts of various cell lines, using IL1B antibody .Secondary antibody: HRP Goat Anti-Rabbit IgG at 1:10000 dilution.Lysates/proteins: 25ug per lane.Blocking buffer: 3% nonfat dry milk in TBST.



Western blot - IL1B Polyclonal Antibody. Western blot analysis of extracts of various cell lines, using IL1B antibody at 1:1000 dilution.Secondary antibody: HRP Goat Anti-Rabbit IgG at 1:10000 dilution.Lysates/proteins: 25ug per lane.Blocking buffer: 3% nonfat dry milk in TBST.Exposure time: 10s.