

Anti-Rel B Mouse antibody

Catalog Number: 169128

Product name

Anti-Rel B Mouse antibody

Specificity

Human

Antibody description

Rel B Mouse monoclonal antibody

Preparation

Antigen: Purified recombinant fragment of human RELB (AA: 273-443) expressed in E. Coli.

Formulation

Purified antibody in PBS with 0.05% sodium azide.

Storage

Store at 4°C short term. Store at -20°C long term.
Avoid freeze / thaw cycle.

Clonality

Monoclonal


Ig Type

IgG2a


Applications

IHC-P, ICC/IF, FC


Validations



Immunohistochemical analysis of paraffin-embedded rectum cancer tissues using RELB mouse mAb with DAB staining.



Immunohistochemical analysis of paraffin-embedded cervical cancer tissues using RELB mouse mAb with DAB staining.




mouse mAb with DAB staining.

Western blot analysis using RELB mAb against

Anti-Rel B Mouse antibody




Catalog Number: 169128



HEK293 (1) and RELB (AA)

Western blot analysis using RELB mAb against

Immunofluorescence analysis of MCF-7 cells using RELB mouse mAb (green). Blue



Flow cytometric analysis of MCF-7 cells using RELB mouse mAb (green) and negative control (red).

human RELB (AA)

