

# Anti-LTA antibody

Catalog Number: 175007

---

**Product name**

Polyclonal

Anti-LTA antibody

**Ig Type****Specificity**

IgG

Human, Rat

**Applications****Antibody description**

WB

Rabbit polyclonal antibody to LTA

**Dilutions****Preparation**

WB: 1:1,000

This antigen of this antibody was synthetic peptide within human lt-alpha aa 51-104.

**Validations****Formulation**

Liquid, 1\*PBS (pH7.4), 0.2% BSA, 25% Glycerol.  
Preservative: 0.05% Sodium Azide.

**Storage**

Store at +4°C after thawing. Aliquot store at -20°C or -80°C. Avoid repeated freeze / thaw cycles.

Fig1: Western blot analysis on LT-alpha (TNF- $\beta$ ) using anti- LT-alpha polyclonal antibody

**Clonality**