

Anti-LTA antibody

Catalog Number: 175008

Product name

Polyclonal

Anti-LTA antibody

Ig Type**Specificity**

IgG

Human, Rat

Applications**Antibody description**

WB

Rabbit polyclonal antibody to LTA

Dilutions**Preparation**

WB: 1:1,000

This antigen of this antibody was synthetic peptide within c-terminal human Lt-alpha(c-term).

Validations**Formulation**

Liquid, 1*PBS (pH7.4), 0.2% BSA, 25% Glycerol.
Preservative: 0.05% Sodium Azide.

**Storage**

Store at +4°C after thawing. Aliquot store at -20°C or -80°C. Avoid repeated freeze / thaw cycles.

Fig1: Western blot analysis on LT-alpha(TNF-β) using anti- LT-alpha polyclonal antibody

Clonality