

# Anti-LTA antibody



Catalog Number: 175008

## Product name

Polyclonal

Anti-LTA antibody

## Ig Type

## Specificity

IgG

Human, Rat

## Applications

## Antibody description

WB

Rabbit polyclonal antibody to LTA

## Dilutions

## Preparation

WB: 1:1,000

This antigen of this antibody was synthetic peptide within c-terminal human Lt-alpha(c-term).

## Validations

## Formulation

Liquid, 1\*PBS (pH7.4), 0.2% BSA, 25% Glycerol.  
Preservative: 0.05% Sodium Azide.

## Storage

Store at +4°C after thawing. Aliquot store at -20°C or -80°C. Avoid repeated freeze / thaw cycles.

## Clonality



Fig1: Western blot analysis on LT-alpha(TNF- $\beta$ ) using anti- LT-alpha polyclonal antibody