

# Anti-PROS1 antibody

Catalog Number: 175010

## Product name

Anti-PROS1 antibody

## Specificity

Human

## Antibody description

Rabbit polyclonal antibody to PROS1

## Preparation

This antigen of this antibody was synthetic peptide within c-terminal human vitamin k-dependent protein s.

## Formulation

Liquid, 1\*PBS (pH7.4), 0.2% BSA, 40% Glycerol.  
Preservative: 0.05% Sodium Azide.

## Storage

Store at +4°C after thawing. Aliquot store at -20°C or -80°C. Avoid repeated freeze / thaw cycles.

## Clonality

Polyclonal

## Ig Type

Rabbit IgG

## Applications

WB

## Dilutions

WB: 1:1,000

## Validations




Fig1: Western blot analysis on human serum using anti- Vitamin K-dependent protein S polyclonal antibody.