

# Anti-TMEM132B antibody



Catalog Number: 175013

## Product name

Anti-TMEM132B antibody

## Specificity

Human, Mouse

## Antibody description

Rabbit polyclonal antibody to TMEM132B

## Preparation

This antigen of this antibody was synthetic peptide within n-terminal human tmem132b.

## Formulation

Liquid, 1\*PBS (pH7.4), 0.2% BSA, 40% Glycerol.  
Preservative: 0.05% Sodium Azide.

## Storage

Store at +4°C after thawing. Aliquot store at -20°C. Avoid repeated freeze / thaw cycles.

## Clonality

Polyclonal

## Ig Type

Rabbit IgG

## Applications

WB

## Dilutions

WB:1:500-1:1,000

## Validations




Fig1: Western blot analysis on D3 using anti-TMEM132B polyclonal antibody.