

Anti-TMEM132E antibody



Catalog Number: 175016

Product name

Anti-TMEM132E antibody

Specificity

Human, Mouse

Antibody description

Rabbit polyclonal antibody to TMEM132E

Preparation

This antigen of this antibody was synthetic peptide within n-terminal human tmem132e.

Formulation

Liquid, 1*PBS (pH7.4), 0.2% BSA, 40% Glycerol.
Preservative: 0.05% Sodium Azide.

Storage

Store at +4°C after thawing. Aliquot store at -20°C. Avoid repeated freeze / thaw cycles.

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications

WB

Dilutions

WB:1:500

Validations




Fig1: Western blot analysis on D3 using anti-TMEM132E polyclonal antibody.