

# Anti-FAT4 antibody



Catalog Number: 175019

## Product name

Anti-FAT4 antibody

## Specificity

Human

## Antibody description

Rabbit polyclonal antibody to FAT4

## Preparation

This antigen of this antibody was synthetic peptide within human fat tumor suppressor homolog 4 aa 2260-2313.

## Formulation

Liquid, 1\*PBS (pH7.4), 0.2% BSA, 40% Glycerol. Preservative: 0.05% Sodium Azide.

## Storage

Store at +4°C after thawing. Aliquot store at -20°C or -80°C. Avoid repeated freeze / thaw cycles.

## Clonality

Polyclonal

## Ig Type

Rabbit IgG

## Applications

WB, ICC

## Dilutions

WB:1:500-1:1,000

ICC: 1:50-1:100

## Validations




Fig1: Western blot analysis on recombinant FAT4 using anti-FAT4 polyclonal antibody.




Fig2: ICC image of FAT4 antibody stained D3 (Murine Embryonic Stem Cell Line).The secondary antibody (red) was goat anti-rabbit IgG (H+L) alexa555 conjugated.