

# Anti-Ngfr antibody

Catalog Number: 175024

## Product name

Anti-Ngfr antibody

## Specificity

Mouse, Rat

## Antibody description

Rabbit polyclonal antibody to Ngfr

## Preparation

This antigen of this antibody was peptide

## Formulation

Liquid, 1\*PBS (pH7.4), 0.2% BSA, 40% Glycerol.  
Preservative: 0.05% Sodium Azide.

## Storage

Store at +4°C after thawing. Aliquot store at -20°C or -80°C. Avoid repeated freeze / thaw cycles.

## Clonality

Polyclonal

## Ig Type

IgG

## Applications

WB, ICC, IHC-P, FC

## Dilutions

WB: 1:2,000~1:5,000

## Validations




Fig1: Western blot analysis on D3 (A) and PC12 (B) using anti-CD271 polyclonal antibody.