

Anti-Lnpk antibody

Catalog Number: 175031

Product name

Anti-Lnpk antibody

Specificity

Mouse, Rat

Antibody description

Mouse monoclonal antibody to Lnpk

Preparation

This antigen of this antibody was recombinant protein

Formulation

Liquid, 1*PBS (pH7.4), 0.2% BSA, 50% Glycerol.
Preservative: 0.05% Sodium Azide.

Storage

Store at +4°C after thawing. Aliquot store at -20°C or -80°C. Avoid repeated freeze / thaw cycles.

Clonality

Monoclonal

Ig Type

Mouse IgG1

Applications

WB, IHC-P

Dilutions

WB: 1:500

IHC-P: 1:50-1:100

Validations





Fig1: Western blot analysis of Lnpk on mouse brain tissue lysate using anti-Lnpk antibody at



1/100 dilution.

Fig2: Immunohistochemical analysis of paraffin-embedded rat spinal cord tissue using anti-Lnpk antibody. Counter stained with hematoxylin.




Fig3: Immunohistochemical analysis of paraffin-embedded mouse brain tissue using anti-Lnpk antibody. Counter stained with hematoxylin.