

# Anti-IGLC2 antibody

Catalog Number: 175051

## Product name

Anti-IGLC2 antibody

## Specificity

Human

## Antibody description

Mouse monoclonal antibody to IGLC2

## Preparation

This antigen of this antibody was recombinant protein

## Formulation

Liquid, 1\**TBS* (pH7.4), 1%*BSA*, 40%*Glycerol*.  
Preservative: 0.05% *Sodium Azide*.

## Storage

Store at +4°C after thawing. Aliquot store at -20°C or -80°C. Avoid repeated freeze / thaw cycles.

## Clonality

Monoclonal

## Ig Type

Mouse IgG1

## Applications

WB, IHC-P

## Dilutions

WB: 1:500

IHC-P: 1:50-1:200

## Validations






Fig1: Western blot analysis of IGLC2 on human IGLC2 recombinant protein using anti-IGLC2



antibody at 1/1,000 dilution.

Fig2: Western blot analysis of IGLC2 on HEK293 (1) and IGLC2-hIgGFc transfected HEK293 (2) cell lysate using anti-IGLC2 antibody at 1/1,000



dilution.

Fig3: Immunohistochemical analysis of paraffin-embedded human colon tissues using anti-IGLC2 antibody. Counter stained with hematoxylin.




Fig4: Immunohistochemical analysis of paraffin-embedded human endometrial cancer tissue using anti-IGLC2 antibody. Counter stained with hematoxylin.