

Anti-FCSK antibody



Catalog Number: 175063

Product name

Anti-FCSK antibody

WB: 1:500-1:1,000

FC: 1:100-1:200

Specificity

Human, Mouse, Rat, Monkey

Antibody description

Mouse monoclonal antibody to FCSK

Preparation

This antigen of this antibody was recombinant protein

Formulation

Liquid, 1*TBS (pH7.4), 1%BSA, Preservative: 0.05% Sodium Azide.

Storage

Store at +4°C after thawing. Aliquot store at -20°C or -80°C. Avoid repeated freeze / thaw cycles.

Clonality

Monoclonal

Ig Type

Mouse IgG1

Applications

WB, FC

Dilutions

Validations




Fig1: Western blot analysis of Fucokinase on different cell lysates using anti-Fucokinase antibody at 1/1,000 dilution.; Positive control; Lane 1: Hela; Lane 2: HepG2; Lane 3: Jurkat; Lane 4: A431; Lane 5: HEK293; Lane 6: MCF-7; Lane 7: PC-12; Lane 8: Cos7; Lane 9:NIH/3T3




Fig2: Flow cytometric analysis of Hela cells with Fucokinase antibody at 1/100 dilution (green) compared with an unlabelled control (cells without incubation with primary antibody; purple).