

Anti-MSX1 antibody



Catalog Number: 175065

Product name

Anti-MSX1 antibody

Specificity

Human

Antibody description

Mouse monoclonal antibody to MSX1

Preparation

This antigen of this antibody was recombinant protein

Formulation

Liquid, 1*TBS (pH7.4), 1%BSA, Preservative: 0.05% Sodium Azide.

Storage

Store at +4°C after thawing. Aliquot store at -20°C or -80°C. Avoid repeated freeze / thaw cycles.

Clonality

Monoclonal

Ig Type

Mouse IgG1

Applications

WB

Dilutions

WB: 1:500-1:1,000

Validations




Fig1: Western blot analysis of MSX1 on NTERA-2 cell lysate using anti-MSX1 antibody at 1/1,000 dilution.