

Anti-H1-1 antibody

Catalog Number: 175100

Product name

Anti-H1-1 antibody

Specificity

Human, Mouse, Rat

Antibody description

Mouse monoclonal antibody to H1-1

Preparation

This antigen of this antibody was synthetic peptide of histone h1 (mono methyl lys25)

Formulation

Liquid, PBS, pH 7.4, containing 0.02% sodium azide as Preservative and 50% Glycerol.

Storage

Store at +4°C after thawing. Aliquot store at -20°C or -80°C. Avoid repeated freeze / thaw cycles.

Clonality

Monoclonal

Ig Type

Mouse IgG

Applications

WB

Dilutions

WB: 1:500-1:1000

Validations




Fig1.; Western blot analysis of 1) HeLa, 2) 3T3, diluted at 1:2000.. Secondary antibody (catalog#: HA1001) was diluted at 1:20000